



# Coag2D

**THE MOST VERSATILE 2 CHANNEL  
OPTICAL COAGULOMETER  
AT YOUR HAND**

**small compact design  
supports D-dimer and AT-III test  
ready for coagulation, chromogenic and turbidimetric tests**



# Coag 2D

## COAG 2D COAGULOMETER SPECIFICATIONS

### Key features:

- 2 independent measuring channels
- 2 reagent positions (2 mL, 5mL)
- 4 incubation positions
- Coagulation screening test
- Quantitative D-dimer measurement
- AT-III activity measurement
- QC feature
- Supports application of Master Calibration Curve
- Maintenance free
- Up to 1000 results storage
- External data archiving possibility
- Connection to Coag Data data management software

### System elements:

- 2 tempered reagent positions with magnetic stirring for heating and mixing
- Colour LCD touchscreen
- Innovative user interface
- LIS interface connection
- Small, compact footprint, easy too operate
- 405 nm, 570 nm and 640 nm wavelength detection
- Optional BCR and printer



#### EUROPE

Diagon Ltd - Hungary  
1047 Budapest, Baross u. 48-52.  
phone.: (+36 1) 369 6500  
fax: (+36 1) 369 6301  
e-mail: [diagon@diagon.com](mailto:diagon@diagon.com)

Diagon GmbH – Austria  
GZO – Zukunftsweg 1.  
Siegendorf, Burgenland 7011  
phone: (+43 2687) 42090  
e-mail: [diagongmbh@diagon.com](mailto:diagongmbh@diagon.com)  
[www.diagon.com](http://www.diagon.com)

#### SOUTH AMERICA

Diagon Ltda - Brasil  
Av. Joaquim José Diniz, 787  
CEP 31910-520  
Bairro Fernão Dias  
Belo Horizonte/MG  
phone: (+55 31) 3488 6777  
fax: (+55 31) 3488 6778  
e-mail: [diagon@diagon.com.br](mailto:diagon@diagon.com.br)  
[www.diagon.com.br](http://www.diagon.com.br)

#### ASIA

PT Inti Diagonama Selaras - Indonesia  
Perkantoran Plaza Pasifik A3/59 -  
Jl. Raya Boulevard Barat  
Kelapa Gading Permal Jakarta 14250  
phone.: (+62 21) 458 40649, 458 40650  
fax: (+62 21) 458 45024  
e-mail: [dsu@centrin.net.id](mailto:dsu@centrin.net.id)