



Coag4D

4-CHANNEL COAGULOMETER

visual sample traffic guidance
supports D-dimer and AT-III tests
ready for coagulation, chromogenic and turbidimetric tests



Coag 4D

4-CHANNEL COAGULOMETER

SEMI-AUTOMATED COAGULATION INSTRUMENT

Key features:

- 4 independent measuring channels (2 parallels or 4 single tests)
- Full-scope coagulation testing
- (Screening tests, Factors, Inhibitors & User defined tests)
- Quantitative D-dimer measurement
- AT-III activity measurement
- Maintenance free
- Liquid crystal color touch screen
- Innovative user interface
- Visual sample traffic guidance for easy and touchless operation
- Synchronized on screen and LED real-time channel status follow-up
- QC feature
- On-board storage of up to 20 cuvettes for preincubation

- Green for preparation
- Orange for incubation
- Red for measurement



System elements:

- 2 thermostated reagent positions both with magnetic stirring
- Easy to use thermal printer
- 2 USB connections for software upgrade and data transfer
- Optional barcode reader
- 405 nm, 570 nm and 640 nm wavelength detection



EUROPE

Diagon Ltd - Hungary
1047 Budapest, Baross u. 48-52.
phone.: (+36 1) 369 6500
fax: (+36 1) 369 6301
e-mail: diagon@diagon.com

Diagon GmbH – Austria
GZO – Zukunftsweg 1.
Siegendorf, Burgenland 7011
phone: (+43 2687) 42090
e-mail: diagonmbh@diagon.com
www.diagon.com

SOUTH AMERICA

Diagon Ltda - Brasil
Av. Joaquim José Diniz, 787
CEP 31910-520
Bairro Fernão Dias
Belo Horizonte/MG
phone: (+55 31) 3488 6777
fax: (+55 31) 3488 6778
e-mail: diagon@diagon.com.br
www.diagon.com.br

ASIA

PT Inti Diagonama Selaras - Indonesia
Perkantoran Plaza Pasifik A3/59 -
Jl. Raya Boulevard Barat
Kelapa Gading Permal Jakarta 14250
phone.: (+62 21) 458 40649, 458 40650
fax: (+62 21) 458 45024
e-mail: dsu@centrin.net.id

